

## Human Lung Directed Differentiation Protocol

F. Hawkins

7/17/15

### 1) Definitive endoderm induction (Day 0 to Day 3-4):

- A. Generate definitive endoderm following the stem cell technologies STEMdiff Definitive Endoderm kit protocol (see attached) in a 6 well plate. Check C-kit/CXCR4 co-expression to assess sufficient definitive endoderm induction. This needs to be assessed for each iPSC/ESC and can change over time. Typically the majority of cells are C-kit+/CXCR4+ between 72 and 96 hours.

### 2) Anterior foregut endoderm (typically 72 hours):

- A. Prepare 12ml of "DS/SB" media per well of definitive endoderm. Divide equally between two conicals and add Y-27632 (ROCK inhibitor) to one (final concentration of 10uM). Place "DS/SB + Y-27632" at room temp. Label and place the other conical of "DS/SB" at 4 degrees celcius. Prepare fresh matrigel coated plate (see Stem Cell technologies E8 or mTeSR1 protocol).
- B. Aspirate definitive endoderm media and gently wash with 1ml of CSFDM. Add 1ml/well of Gentle Cell Dissociation Reagent. Wait 2-3 min. Aspirate. Add 1ml "DS/SB + Y-27632" directly to the well. Gently triturate. The cells should very easily detach from the tissue culture plate as clumps. Plate the cells at a density of 200K cells/cm<sup>2</sup> in "DS/SB + Y-27632" on freshly coated matrigel plates.
- C. The following day change the media to "DS/SB". The duration of "DS/SB" may need to be optimized for NKX2-1 induction. In our experience, 48 hours is insufficient. For most iPSC/ESC we use 72 hours.

### 3) NKX2-1 Induction (Day 6 to 15):

- A) Prepare "CFKBRa" media with freshly prepared retinoic acid.
- B) Aspirate "DS/SB" media.
- C) Add 2ml/well of "CFKBRa".
- D) Refeed every 24-48 hours.

### 4) Expansion (Day 15-22):

- A) Re-plating unsorted cells on matrigel:
  - I. Prepare fresh matrigel coated plates (see Stem Cell technologies E8 or mTeSR1 protocol).
  - II. Prepare "CFK media"
  - III. Aspirate "CFKBRa" media and wash X2 with CSFDM.
  - IV. Add warm 0.05% trypsin and place the tissue culture plate in the incubator for 2-4 minutes.
  - V. Aspirate trypsin and wash X 1 with STOP media (DMEM +10% FBS) and X1 with CSFDM.

- VI. Add 1ml of "CFK" and triturate gently until the cell sheet has been dissociated into small clumps.
- VII. Transfer to a 15 ml conical. Dilute to a 1:3 split ratio i.e. 1 well of a 6 well plate should be re-suspended with 6 ml of "CFK" and supplemented with Y-27632 (final concentration 10uM) and divided between three wells of freshly coated 6 well plate.
- VIII. Replace media every 24-48 hours.

B) Re-plating unsorted cells *in* 3-D matrigel drops:

- I. Thaw matrigel (356231) on ice
- II. Repeat steps I to V above.
- III. Add 1ml of CSFDM/well and triturate gently until the cell sheet has been dissociated into small clumps.
- IV. Divide the dissociated 1ml sample between 1.5ml eppendorfs (approx. 50ul per eppendorf). Add 1ml CSFDM per eppendorf and centrifuge at 200G for 5 min.
- V. Aspirate as much of the supernatant as possible without displacing the cell pellet (to avoid diluting the matrigel).
- VI. Place eppendorf on ice for 2 min.
- VII. Re-suspend each cell pellet in 50ul of matrigel using cold 200ul pipette tips. Take care to (a) avoid bubbles and (b) distribute the clumps evenly in the matrigel.
- VIII. Slowly pipette the matrigel + cells into the center of a 12 well plate (50uL/drop).
- IX. Place 12 well plate in the incubator for 20 min.
- X. Add CFK + 10uM Y-27632 media.
- XI. Re-feed with CFK every 24-48 hours.

C) Re-plating sorted NXX2-1+ cells in 3-D matrigel drops:

- I. Prepare FACS buffer:
  - Hank's Balanced Salt Solution
  - 2% FBS
  - Primocin
  - HEPES 25mM
  - EDTA 2mM
- I. Wash day 15 cells with CSFDM
- II. Use a 10ul pipette to etch many scrapes across the monolayer
- III. Add 1ml warm 0.05% trypsin and place in incubator for 14-18 minutes
- IV. Triturate (gently) until cell sheet is detached from plate and transfer to a 15ml conical
- V. Add an additional 1ml of trypsin. Manually shake, flick and roll the 15ml conical for 3-5 minutes. Stop occasionally and allow clumps to settle for 20-30 sec. Transfer cloudy supernatant to a separate 15ml conical containing "STOP Media". Continue shaking, flicking and rolling the remaining clumps until mostly dissociated and then transfer the remaining cells to the conical containing "STOP Media".
- VI. Centrifuge at 200G x 5min.

- VII. Aspirate supernatant and re-suspend in FACS buffer
- VIII. Filter through a 40um filter x 2.
- IX. Add Y-27632 (final concentration 10uM) and Propidium Iodide to FACS sample and prepare collection tubes with FACS buffer supplemented with Y-27632 (final concentration 10uM).
- X. Sort NKX2-1+/PI negative
- XI. Spin collected samples at 300G for 7min
- XII. Repeat steps B(V)-B(XI) above.

#### 5) Maturation (Day 22+):

- I. On day 22 prepare "CFK+DCI" media
- II. Change media every 24-48 hours