

Hemoglobin Variant Analysis Using MALDI-reISD

Roger Th  berge¹, Sergei Dikler², Mark E. McComb¹, Catherine E. Costello¹

¹Center for Biomedical Mass Spectrometry, Boston University School of Medicine, Boston, MA. ²Bruker Daltonics

Novel Aspect

Top-Down analysis of hemoglobin variants using MALDI-ISD

Introduction – Limit 120 words

More than 1200 recognized hemoglobin structural variants exist in the human population, and many of these underlie phenotypic diseases. Mass spectrometry constitutes a rapid and accurate means for the detection and characterization disease-specific protein variants. In our continuing efforts to develop mass spectrometry based methods to detect and characterize hemoglobin variants, we report preliminary results using MALDI-reISD (in-source decay) for the purpose of hemoglobin variant analysis.

Method – Limit 120 words

Whole blood was diluted approximately 1:250 in water. The sDHB matrix solution was prepared to a concentration of 50 g/L in 50% acetonitrile/water/0.1 formic acid. A saturated solution of 1,5-diaminonaphthalene (1,5-DAN) matrix was prepared in 50% acetonitrile/ water/0.1 formic acid. Analyses were carried out for target spots that were dried down from equal volumes of sample and matrix solutions, using 1 μ L of each. MALDI-reISD spectra were acquired on a Bruker UltrafleXtreme MALDI-TOF/TOF MS and consisted of signals generated from 12000 to 15000 accumulated laser shots. The resulting fragment masses were analyzed using BioTools and BUPID-Topdown (Boston University Protein Identifier-Topdown), a custom-programmed software algorithm written in-house. Caution: 1,5-DAN is a carcinogen.

Preliminary data – Limit 300 words (299)

The MALDI-reISD spectrum of normal hemoglobin in sDHB exhibited extensive fragmentation consisting mainly of c, z and y-ions describing the N- and C-termini of both the alpha and beta chains. Although both y and z-ion series were present, the signals from the y-ion series were more intense and extensive. The fragmentation roughly covered the first 50 amino acids of the N- and C-termini. In the analysis of a sickle cell sample, a -30 Da peak series was observed for the N-terminal fragmentation of the beta chain, consistent with the presence of the Glu6Val mutation. An alpha chain variant, hemoglobin Westmead (His122Gln), was also analyzed and characterized using MALDI-reISD. We have also obtained data using 1,5-DAN as the matrix and have noted a few significant differences between the spectra recorded with 1,5-DAN and the sDHB data. Generally, 1,5-DAN favors N-terminal fragmentation and generates higher

fragment ion abundances than does sDHB. Direct analyses of hemoglobins from blood samples are well suited for MALDI-LS/MS since sample purity and quantity are not issues. The "purification" of Hb samples can be performed by simple dilution of whole blood (Hb conc. 124 g/L), whereby the abundances of other blood components become almost insignificant with respect to hemoglobin. Also, sample amount is not an issue and this factor is generally cited as a significant hurdle to applications of MALDI-LS/MS in proteomics. The main limitation of the technique at this point is the lack of coverage of the core region of the alpha and beta chain sequences. MALDI-reLS/MS offers a simple, fast and potentially high throughput technique for the detection of hemoglobin variants. We are currently working to optimize analysis conditions..

This research is supported by NIH P41 RR10888. The authors are grateful to the BUSM sickle cell center for providing the blood samples with ISB approval.